

FRAP Analysis Report

Fluorescence Recovery After Photobleaching

Samples: 100 μ M OA + SCDi vs SCDi

Date: December 12, 2025

Executive Summary

This report presents FRAP analysis comparing two treatment conditions: 100 uM Oleic Acid + SCDi and SCDi alone.

Key Findings

- 1. Recovery Kinetics: SCDi alone shows slower recovery ($t_{1/2} = 1.31s$) compared to 100 uM OA + SCDi ($t_{1/2} = 0.60s$), representing a ~2.2x difference.
- 2. Mobile Fractions: Both treatments show similar mobile fractions (~55%).
- 3. Model Fit: Single exponential provides excellent fits ($R^2 > 0.99$).

Summary Statistics

Parameter	100 uM OA + SCDi	SCDi	Difference
Mobile Fraction (%)	55.6	54.5	-1.2
Half-time $t_{1/2}$ (s)	0.601	1.306	0.705
Time Constant τ (s)	0.867	1.884	1.017
R-squared	0.9998	0.9962	-

Figures

Figure 1: Main Publication Figure





Figure 1. FRAP analysis of 100 μ M OA + SCDi and SCDi treatments. (A) Representative images. (B) Recovery curves. (C) Mobile fractions. (D) Half-times. (E) Summary statistics.

Figure S1: Individual Recovery Curves



Conclusions

This FRAP analysis reveals distinct molecular dynamics between treatments:

1. SCDi alone results in slower fluorescence recovery ($t_{1/2} = 1.31\text{s}$) compared to 100 μM OA + SCDi ($t_{1/2} = 0.60\text{s}$).
2. Both treatments show similar mobile fractions (~55%).
3. The faster recovery with OA + SCDi suggests oleic acid may increase molecular mobility despite SCD inhibition.